



Cuticell®



Protecting,
soothing,
helping to heal:
our **ointment**
dressing

The ointment dressing with Eucerit®.




Cuticell® is soothing and low-adherent, thus minimizing pain at dressing change. The primary wound contact layer is made of smooth acetate fabric impregnated with non-medicated ointment on Eucerit® basis. The result: an ointment dressing that ensures that the wound is not disturbed and helps maintain a moist environment that supports wound healing.

Special characteristics



Impregnation with ointment on Eucerit® basis

Helps maintain a moist environment to support wound healing. Does not charge the wound with excess ointment.



Smooth acetate fabric with knit structure

Lint-free and does not fray when cut to size. The defined small mesh size of 1 x 1 mm prevents tissue ingrowth.



Permeable to air and exudate

Allows exudate to drain freely into an absorbent secondary dressing. Prevents occlusive effects.

Advantages

- Low-adherent and soothing
- Minimizes pain at dressing change
- Supports a moist environment
- Does not disturb wound healing
- Also suitable for patients with sensitive skin
- Simple handling without release paper

Indications

For superficial exuding wounds such as

- First and second-degree burns
- Cuts and abrasions
- Lacerations
- Radiation injuries
- Donor and recipient skin graft sites

Range available

sterile, individually sealed

REF No.	Size	Pack size	Shipping Unit
72539-00	7.5 cm x 7.5 cm / 3 in. x 3 in.	10	600
72539-01	7.5 cm x 7.5 cm / 3 in. x 3 in.	50	600
72539-02	7.5 cm x 20 cm / 3 in. x 8 in.	10	300
72539-03	7.5 cm x 20 cm / 3 in. x 8 in.	50	300
72539-04	20 cm x 40 cm / 8 in. x 16 in.	25	150

