

WEREWOLF* COBLATION* System introduces the new innovative FLOW 90 wand



The WEREWOLF* COBLATION* System can be used on a range of joint tissue types to provide COBLATION technology across a range of different settings.

Designed for surgeons who want to restore joints as close as possible to their natural anatomy.





Subacromial decompression





Labral Peeling

Technical specifications

Controller	WEREWOLF Controller
Shaft diameter	3.76 mm
Working length	137 mm
Tip width	3.40 mm
Default set point	Hi (•)

Ordering Information

WEREWOLF COBLATION FLOW 90° WAND	
Reference #	Description
72290105	WEREWOLF Generator
72290038	FLOW 90 wand

For more information, contact your sales representative or visit smith-nephew.com.