


## Cuticell® Classic



Protecting,  
soothing,  
helping to heal:  
our **paraffin gauze**

# Cuticell® Classic

## The paraffin-impregnated dressing.




Cuticell® Classic is a low-adherent wound contact layer made of open weave cotton gauze impregnated with soft paraffin – to facilitate dressing change and to support a moist environment for the wound.

### Special characteristics



#### Impregnation with soft paraffin

Makes for low adherence and helps maintain a moist environment to support wound healing.



#### Interlocking weave gauze fabric

For increased stability when cut to size.  
Does not fray.



#### Permeable to air and exudate

Allows exudate to drain freely into an absorbent secondary dressing.  
Prevents occlusive effects.

### Advantages

- Low-adherent and soothing
- Facilitates dressing change
- Supports a moist environment
- Permeable to air and exudate
- Also available in handy dispenser box for economical treatment

### Indications

For exuding wounds such as

- First and second-degree burns
- Lacerations and abrasions
- Leg ulcers
- Donor and recipient skin graft sites

### Range available

sterile, individually sealed

REF No.	Size	Pack size	Shipping Unit
72538-00	5 cm x 5 cm	5	600
72538-01	5 cm x 5 cm	50	600
72538-02	10 cm x 10 cm	10	600
72538-03	10 cm x 10 cm	100	600
72538-04	10 cm x 40 cm	10	300

sterile, plastic tin or dispenser roll

REF No.	Size	Pack size	Shipping Unit
72538-05	10 cm x 10 cm	36	12 tins
72538-06	10 cm x 7 m	1	12 tins
72538-07	15 cm x 2 m	12	24 rolls



70550-00341-00 03/2008